



Biological Safety Cabinet ClassII A2

| Biological Safety Cabinet

Details Info

- Features:
- The motorized glass sash can be lift exactly to the 20cm height. The safety for opening only by 'one press' on the control panel
  - The HEPA filter can be replaced by easy and safe way from the front
  - The Colorized Visual LCD Display indicate the status of the cabinet with a simple and convenient international figures and code Elegance R70 are of the front panel to avoid any hurt by sharp edge.
  - Adopts the spring gas bar mechanism to open the glass door frame and front panel automatically , save operator's time and easy for clean and maintenance choose the unique shape inflow baffle to assure the airflow average and have the same function as the hand rack
  - Compare to it's overall size ,the work area offers a big space for operator .Adopts universal wheels make the cabinet can remove smooth and easily . Two airflow sensors for inflow and down flow mounted interior

Technical Data:

HEPA Filters	Two Hepafilter for intake and exhaust air
Filter Efficiency	The filters with efficiency of 99.9995%at0.3/um
Material	Interior :304 stainless steel .Exterior: cold –rolled steel with epoxy coating
Velocity	Down flow :0.3m/s inflow:0.45m/s
Airflow Mode	30% for exhaust.70%forcirculation
Noise	≤60dB
Display	Colorized LCD display
Light	1000 IUX UV: 30w
Power Supply	220V 50HZ

Dimensions:

Model	BSC-1050IIA2	BSC-1350IIA2	BSC-1650IIA2	BSC-1950IIA2
Overall size(mm)W*D*H)	1050*750*1990	1350*750*1990	1650*750*1990	1950*750*1990
Internal size(mm)W*D*H)	910*590*650	1203*590*650	1500*590*650	1800*590*650

Head Office :1600 AVENUE DE LORIMIER BUREAU, 384 MONTRÉAL (QUÉBEC) CANADA H2K3W8 USAOffice:128 W SUFFLOK AVE,CENTRAL ISLIP,11722 N Middle east & Africa  
Office : warehouse 9 /shed 9, Technology park P.O. Box 54806 RAK FZE , UAE

Email Us At:  
sales@acculabusa.com middleeast@acculabusa.com

Phone Support: